


NANOVEA ST400

***STANDARD
OPTICAL PROFILER***



microworld®

NANOVEA.COM

SETTING STANDARDS

Designed with Chromatic Light technology, which measures physical wavelength, ST400 Profiler provides the highest accuracy on any roughness, any form, any material. Transparent or opaque.

HIGHLY CUSTOMIZABLE

WIDE RANGE OF ADD-ONS

ROTATIONAL OPTIONS


STAGE OR CYLINDER

ULTRA FAST MEASUREMENTS


w/ HIGH SPEED LINE SENSOR

up to 200 mm

HEIGHT SAMPLE CLEARANCE



With larger X-Y stages, 360° rotational stages & many other custom configurations available, ST400 is ideal for the widest range of samples with varied geometries.



X-Y
STAGE TRAVEL
200 x 150 mm

Z
AXIS
50 mm Motorized


X-Y
MAX SPEED
40 mm/s

THE POWER OF CHROMATIC LIGHT


*NANOVEA's Non-Contact Optical Profilers are the ideal upgrade
from traditional contact stylus and laser profilometers.*




2D & 3D NON-CONTACT SURFACE MEASUREMENTS




ROUGHNESS & FINISH




TEXTURE




GEOMETRY & SHAPE



STEP HEIGHT & THICKNESS



VOLUME & AREA



FLATNESS & WARPAGE

1 nm MAX VERTICAL RESOLUTION


up to 87° MAX SURFACE ANGLE



STANDARD SENSOR

SINGLE POINT

	PS1	PS2	PS3	PS4	PS5	PS6
MAX HEIGHT RANGE	110μm	300μm	1.1mm	3.5mm	10mm	24mm
WORKING DISTANCE	3.35mm	10.8mm	12.0mm	16.2mm	25.9mm	20mm
LATERAL X-Y ACCURACY	0.9μm	1.2μm	2.0μm	3.0μm	7.0μm	8.0μm
HEIGHT REPEATABILITY*	1.2nm	2.2nm	3.4nm	17nm	31nm	41nm



HIGH SPEED SENSOR

192 POINTS


	LS1	LS2	LS3
MAX HEIGHT RANGE	200μm	0.95mm	3.9mm
WORKING DISTANCE	5.3mm	18.5mm	41mm
HEIGHT REPEATABILITY Ra*	14nm	21nm	70nm
LINE WIDTH	0.96mm	1.91mm	4.78mm
PITCH	5μm	10μm	25μm
LATERAL ACCURACY OF EACH POINT	1μm	2μm	5μm
ACQUISITION RATE (points per second)	384KHz	384KHz	384KHz

* Fixed point measurement on glass. Ra average height variation for 1,200 points (100 sampling).

HOW IT WORKS

Chromatic Light Technology operates via a process that utilizes white light and a series of spherochromatic lenses. The spherochromatic lenses split the white light into individual wavelengths with unique vertical focal points (vertical distance from surface or height). All wavelengths and their corresponding heights make up the height range measurement scale of a sensor.

The wavelength with the highest intensity will be detected by the spectrometer which processes the wavelength's associated height. During a full raster scan, this process takes a fraction of a second and produces an accurate height map of the surface of interest.



NO COMPLEX ALGORITHMS

NO DATA STITCHING WITHIN X-Y STAGE TRAVEL

THE PROBLEM WITH OTHER TECHNIQUES

LATERAL RESOLUTION vs LATERAL ACCURACY

1 px



3 x 3 px



6 x 6 px



30 x 30 px



135 x 135 px



520 x 520 px



NOT ENOUGH DATA TO CALCULATE FOCUS
NO PRACTICAL USE

PIXEL SIZE RESOLUTION: 2nm

FOCUS CAN BE CALCULATED
THE SMALLEST INCREMENT
FOR ANY PRACTICAL USE

EFFECTIVE ACCURACY: 1040 nm

THEM

Camera Pixel Size or **Display Resolution** is often defined as **lateral resolution** to impress clients. Instruments that use camera pixel-based technology require complex algorithms to determine the focal point of the instrument which is problematic for complex surfaces.

US

Chromatic Light provides **lateral accuracy** which is determined by physics and is directly related to the spot size of the chromatic light source of the optical sensor.

LASER SCANNING CONFOCAL MICROSCOPE



LASER RADIATION

HEALTH HAZARD

Exposure to laser light reflectivity

INCONSISTENT LASER LIGHT WAVELENGTH

Inconsistencies in wavelength during scanning
affect accuracy of results

DECEPTIVE 'DISPLAY RESOLUTION'

Lateral & height accuracy are fixed by the objective lens
making 'Display Resolution' insignificant

COMPLEX ALGORITHMS

Alpha blending algorithms stitch collected data
layer by layer, grounding accuracy on complex calculations

STITCHING REQUIRED

Objective lenses have limited fixed fields of view
Stitching of larger areas compromises accuracy of the scan

50x SLOWER

Data acquisition speed up to 7.9 KHz

VS

CHROMATIC LIGHT OPTICAL SENSOR

SAFE WHITE LIGHT

No need for protective wear

UNIFORM & BROAD WHITE LIGHT SPECTRUM

Changes in wavelength are the data being collected

INDEPENDENT LATERAL & HEIGHT ACCURACY

Lateral & height accuracy can be mixed and matched
to meet a broad range of scanning requirements

NO ALGORITHMS

Physical wavelength reflected from the surface
is measured directly for an accurate representative height map

NO STITCHING

Data points are collected continuously providing
the same level of accuracy for both small and large areas

50x FASTER

Data acquisition speed up to 384 KHz

LATERAL ACCURACY

For 50x objective (370 x 277 μm)

$\pm 2\%$ of measuring value

$\pm 2\% \times 370 \mu\text{m}$

$\approx 15 \mu\text{m}$

w/ stitching algorithms $\gg 15 \mu\text{m}$



Step size:

= $5 \mu\text{m}$

3x BETTER LATERAL ACCURACY

HEIGHT ACCURACY

$\approx 0.2 + L/100 \mu\text{m}$

$\approx 0.2 + 950/100 \mu\text{m}$

$\approx 9.7 \mu\text{m}$



950 μm range

$\approx 0.6 \mu\text{m}$

16x BETTER HEIGHT ACCURACY

AREA TESTED

STITCHING REQUIRED

scans (25 x 25 mm)

25 000 μm / 370 μm x 25 000 μm / 277 μm

68 x 91

= **6188 scans**



NO STITCHING

Consistent accuracy across any measurement size

1 SCAN


TEST TIME

6 sec per scan

+ 4 sec displacement & stitching

= 10 sec/scan x 6188 scans

= **61860 seconds** (≈ 17 hours)




Scan time (25 x 25 mm)

= **29.6 seconds**

2090x FASTER

MICROSCOPE VIDEO IMAGING



up to 100x
objective magnification

1200 x 1600
color video camera

ACCURACY OF <math><0.2 \mu\text{m}</math>
to/from indenter position

LARGE AREA STITCHING
capability

NANOVEA

ST400

OPTICAL PROFILER



microworld®

GRENOBLE - FRANCE

Tel : +33 (0)4 76 56 16 17


Email : contact@microworld.eu

www.microworld.eu

Also available in other configurations



PORTABLE



COMPACT



HIGH SPEED



LARGE AREA



ZERO NOISE

